lacrifresh® comfort

Soothes ocular
dryness, irritations
and eye strain

COMFORT AND WETTING SOLUTION FOR ALL CONTACT LENSES

USAGE:

- **Avizor Lacrifresh Comfort** are lubricating drops that act instantly.
- With similar characteristics to the human tear, they soothe ocular tiredness and dryness.
- They are especially recommended for contact lens wearers that spend long hours in front of a computer screen and/or in stuffy atmospheres (long plane or car trips, air-conditioned rooms, smoky environment, etc.).
- They are compatible with all types of contact lenses.

INSTRUCTIONS FOR USE:

1. Wash your hands with soap and water before using the product.
2. Direct use: Apply one or two drops into each eye and blink several times.
3. Indirect use: Apply one or two drops into the internal part of each lens before putting them on.
4. Use **Avizor Lacrifresh Comfort** as many times as needed while wearing your contact lenses.

Last revision date: 2010/04

PRECAUTIONS:

- Do not let any surface come into contact with the nozzle to prevent contamination.
- Do not use the product after it has been opened for over 60 days, nor after the expiry date.
- Do not use the product if the safety seal is missing or broken.
- Keep well closed after each use.
- Do not ingest the solution and keep away from children's reach.
- Keep at room temperature.
- Do not use the product in case you are allergic to any of its components.
- In case of ocular irritation, consult a specialist.
- Before using any ophthalmic medicine, inform your doctor you are a contact lens user.

COMPOSITION:

Isonic, buffered and aqueous solution, Povidone 1%, EDTA 0,1% and Polyhexanide 0,0002%.

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