9464

BAUSCH + LOMB

blink contacts™

Eye Drops

en

Directions:
To soothe dry, irritated, uncomfortable eyes, complete these simple steps with Blink Contacts™ Eye Drops:
STEP 1: Place 1 or 2 drops in affected eye as reactant lenses there is no need to remove them.
For extra comfort with lenses: place 1 or 2 drops on each lens before inserting.
STEP 2: Blink several times to wash the eye. Use as often as required up to 8 times per day.

Contents: A sterile, isotonic, aqueous buffered, preserved solution that contains Sodium chloride (0.005 %), Sodium hyaluronate (0.15 %), Boric acid, Sodium borate, Sodium chloride, Potassium chloride, Calcium chloride, Magnesium chloride, Purified water.

Actions and Clinical Benefits: When used as directed Blink Contacts Eye Drops provides improved ocular comfort for contact lens wearers without affecting quality of vision. Children, pregnant or breastfeeding women were not included in the clinical trial(s) conducted.

Intended Purpose: Soothe dry, irritated, uncomfortable eyes. Provides lubrication and hydration to surface of the eye. Rewet contact lenses.

Indications (Uses): To lubricate and rewet soft (hydrophilic) and rigid gas permeable (RGP) contact lenses, to help relieve dryness, discomfort and irritation that may be associated with lens wear and to cushion lenses by placing a drop on the lens prior to application on the eye.

Intended Patient Population and Users: Users with dry, irritated or uncomfortable eyes with and without contacts.

Contraindications (Reasons Not to Use): Do not use if you are allergic to any ingredient in this product.

Important Safety Information:

- Always follow the directions for use. Failure to follow product directions may lead to vision loss.
- Visit your eye care practitioner regularly.
- Always wash and dry hands before use.
- To avoid contamination or eye injury do not touch bottle tip on any surface and avoid contact with eyes.
- Blink Contacts Eye Drops will not disinfect your lenses.
- Borates (boric acid & sodium borate) are defined as CMR 1B substances in a concentration above 0.1 % weight by weight and are safe when the product is used according to label instructions.
- Consult your doctor or eye care practitioner prior to using in children, pregnant or breastfeeding women.

Precautions:

- Keep the bottle tightly closed when not in use.
- Use before the expiration date.
- Discard any remaining solution 45 days after opening.
- Do not use if imprinted neckband is broken or missing.
- Keep out of reach of children.
- Do not use in conjunction with other eye medications without first consulting your eye care practitioner.
- Dispose as per the local requirement(s).

Storage and Discard Date:

- Store between 5 °C to 25 °C.
- Discard any remaining solution 45 days after opening.

WARNINGS:

- This solution is not intended to disinfect your lenses.
- Using your solution beyond the discard date could result in contamination of the solution and can lead to severe infection, vision loss or blindness.

Adverse Reactions (Problems and What to Do):
You should be aware that the following problems may occur:

- Eyes stinging, burning, itching (irritation), inflammatory reaction or other pain/discomfort
- Abnormal feeling of something in the eye (e.g., foreign body, scratched area)
- Excessive watering (tearing) of the eyes or unusual eye secretions
- Redness of the eyes
- Reduced sharpness of vision (poor visual acuity)
- Sensitivity to light (photophobia)

If you experience these or any other complications, consult with your doctor or eye care practitioner.
If you wear contact lenses and experience these or any other complications, you should immediately remove your lenses. A serious condition such as infection, corneal ulcer, neovascularization or iritis may be present. Seek immediate professional identification of the problem and prompt treatment to avoid serious eye damage.

Reporting of Complaints and Serious Incidents:
Any side effects, adverse reactions, or product complaints from use of this product should be reported to the manufacturer. Country specific information can be found on www.bausch.com/contactus
Reporting of serious incidents from use of this product should be reported to the manufacturer and to the competent authority of your Member State.
Explanation of Symbols located at www.bausch.com/symbols

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