ENGLISH | Instructions for Use

Disinfecting, cleaning and storing system for rigid gas permeable contact lenses.

FOR SCLERAL CONTACT LENSES

Rinse the lenses with **cleadew CareSolution** or other rinsing solutions (sterile saline solution, rinsing solutions for contact lenses, and multipurpose solutions [MPS]) for soft contact lenses) when wearing the lenses after care as required by this lens care system. Avoid using solutions containing sodium chlorite (for example, LENS Plus [Ocupure] or REGARD).

Wash your hands well with mild soap before handling your contact lenses.

Use the special lens case included (do not use any other lens cases).

1 DISINFECTING, CLEANING AND NEUTRALIZING

- Set a neutralizing and cleaning tablet **at the specified position (slit)** of the dedicated case body. Be sure to use the dedicated case. **⚠**Setting the tablet at a position other than the specified position may lead to unsuccessful neutralisation, resulting in colouration of the lens.
- Place the lenses in the lens holder of the case.
- Fill the case with the disinfecting solution up to the marked line.
- Close the case lid and let it stand for at least four hours.
- Confirm that the lenses and the solution are clear.

2 RINSING LENSES

- Pick up the cleaned/disinfected lenses from the holder. Hold them with your fingers and rinse both sides with **cleadew CareSolution** or sterile saline solution for **at least 10 seconds**. **⚠**If you fail to rinse the lenses or your rinsing was not thorough enough, the lenses may cause eye irritation. Be sure to rinse the lenses thoroughly.

STORING LENSES

- To store the lenses, be sure to follow the instructions of eye care professionals.
- When planning to wear lenses that have been stored for 24 hours or longer, it is recommended to start care from the beginning before wearing them.

⚠ PRECAUTIONS

- ⚠**To avoid serious eye infections or damage, always follow the user instructions.
- Do not use the product if you are allergic to any of the ingredients used in it.
- If you experience such symptoms as eye pain and irritation, remove your lenses immediately and consult your eye doctor or an eye care practitioner.
- In the unlikely event that the lens is coloured after care, do not wear the coloured lens.
- Always use the **cleadew SL** system components including the lens case.
- Use the disinfecting solution and the neutralizing and cleaning tablet together.
- Do not use the product together with another solution, tablet, or pharmaceutical drug.
- After using up all of the neutralizing and cleaning tablets, discard any remaining solution in the bottle.
- Do not use for soft contact lenses.
- ⚠**Do not use the product if the bottle or package is soiled or damaged.
- Do not use the tablet if it is damaged or discoloured.
- Do not use the solution if it is cloudy or its colour is a colour other than orange due to discolouration.
- Do not wear lenses without rinsing them with **cleadew CareSolution** or sterile saline solution after disinfecting, cleaning and neutralizing.
- Do not reuse the solution. (The solution after lens care is neutralized and cannot disinfect the lens. To prevent eye infection, do not reuse the solution.)
- Do not use the product after the expiry date.
- ⚠**Discard any remaining solution 90 days after opening the bottle.
- Use the tablets immediately after opening the package.
- Do not touch the tip of the bottle to avoid contamination.
- Close the bottle tightly after use.
- To prevent bacterial contamination, do not transfer the disinfecting solution to another container, nor replenish it with water or other liquid medical agent.
- Do not use tap water to rinse the lenses or lens case.
- ⚠**Keep out of reach of children.
- ⚠**Store products in a dry place at room temperature, avoiding direct sunlight.
- To prevent contamination, use a new lens case with each purchase.
- If the solution adheres to your clothes, wash with water immediately.
- Rinse off any remaining solution on the wash stand with water after each use to prevent staining on the surface of the wash stand.
- Do not swallow the solution or tablet. If you do so, drink a large amount of drinking water immediately.
- ⚠**Not to be used in the eye. If you do so, immediately rinse the eye with sterile saline solution or eye washes.
- Do not change the contact lens care system without consulting your eye care professional.
- ⚠**Boric acid is toxic to reproduction and can reduce future fertility if taken by children or harm babies if taken by pregnant or lactating women. Although this product also contains boric acid, the amount of boric acid in this product is at a concentration that does not pose a health hazard. This product is safe for use by children and pregnant or lactating women as long as the product's directions are followed. Sensitivity to chemicals differs from person to person. For sensitive individuals, the recommendation is to rinse the lens with sterile saline solution before wearing.
- Do not use the solution after the expiry date.
- ⚠**Discard any remaining solution 90 days after opening the bottle.
- Use the tablets immediately after opening the package.
- Do not touch the tip of the bottle to avoid contamination.
- Close the bottle tightly after use.
- To prevent bacterial contamination, do not transfer the disinfecting solution to another container, nor replenish it with water or other liquid medical agent.

INGREDIENTS

[Disinfecting solution]
Povidone-iodine (0.05%)
⚠Boric acid

STERILE A **⚠**

[Neutralizing and cleaning tablet]
Sodium Sulfate (2.4 mg/tablet),
Proteolytic enzyme (8.0 mg/tablet)

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