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Each lens pack has a Lot Number (LOT). If you experience any difficulty with your contact lenses please quote this number.

Daysoft Limited, 5 Livingston Boulevard Hamilton International Technology Park, Blantyre, Scotland G72 0BP

daysoft® is a registered Trade Mark of Daysoft Limited

Advena Limited, Tower Business Centre, 2nd Flr., Tower Street, Swatara, BKR 4013 Malta

Wearer Instructions

Sterile daily disposable contact lenses in a preservative free phosphate buffered saline solution for wearers requiring optical correction of refractive ametropia (myopia and hyperopia) in phakic or aphakic persons with non-diseased eyes, exhibiting astigmatism of 1.00 diopter or less, that does not interfere with visual acuity. Lenses are for use only by wearers that have been assessed as suitable and have been prescribed by an Eye Care Professional (ECP). Ensure you have received instructions in wearing and removing contact lenses by an ECP prior to use.

- Wash and dry your hands before opening the lens pack.
- It is recommended that you see your ECP routinely as directed. Use, insert and remove your contact lenses as in-structed by your ECP. Never insert your contact lenses if you are experiencing eye inflammation, eye dryness or eye disease. If you experience eye discomfort after inserting your contact lenses remove them immediately and contact your ECP as soon as possible.
- If you think that your eyesight has changed, that you require different lens powers or you have any concern about your ocular health, consult your ECP. These contact lenses are not direct replacements for hard (RGP), toric or multifocal lenses. Consult your ECP before changing to or from these lens types.
- These contact lenses are designed for "single-use" daily wear, discard your lenses after removal, do not sleep in your lens, do not attempt to wash or otherwise clean your lens... always use a fresh lens. They should not be reused as they will no longer be sterile, and you risk introducing harmful micro-organisms to the eye.
- Each lens pack has a Lot Number (LOT). If you experience any difficulty with your contact lenses please quote this number.

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Each lens pack also has an Expiry Date (EXP). Never use a contact lens after this date and never use a contact lens if there is any sign of damage to the lens pack.

- Direct exposure of contact lenses to non-sterile water, e.g. tap water, whirlpool bath, swimming, participation in water sports, increases the risk of microbial infection.
- If you use eye-drops when wearing these lenses use those labelled or recommended by your ECP as com-patible with soft contact lenses.
- Some make-up, hairspray and hand creams may cause eye discomfort. To avoid this protect your eyes and con-tact lenses from contact with these products.
- Contact lenses should be removed immediately after contact with noxious vapour, chemical or hazardous substances, or hazardous environments with ocular impact.
- UV absorbing contact lenses are not substitutes for pro-protective UV absorbing eye wear.
- Do not freeze.
- These contact lenses should be worn with a valid prescription and should not be given to friends or family to try without examination by an ECP and valid prescription.
- Should you experience any side effects while using your contact lenses, please discontinue use immediately. Side effects may include bacterial infection, ocular ir-ritation, toxic or allergic reaction, loss of visual acuity and scratches or abrasions on the eye surface.
- Any serious incident that has occurred in relation to these lenses should be reported to the manufacturer and the national competent authority.

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